
SVBONY

(EN:1-4)

(DE:5-8)

(FR:9-12)

(IT:13-16)

(ES:17-20)

(RU:21-24)

(JP:25-28)

SV605

User Manual

EU Importer: Germany Retevis Technology GmbH
Address: Uetzenacker 29,38176 Wendeburg

REFACE

Thank you very much for your purchase of a SVBONY microscope

This manual applies to SV605. Read this manual before using the SV605 microscope.

Setting Up

1. Remove the Styrofoam container from the carton.
2. Carefully remove the microscope and other parts from the container and set them on a table, desk, or other flat surface.
3. Remove eyepieces from the plastic bags.
4. Insert the 10x eyepiece into the eyepiece tube

You are now ready to use your microscope for looking at specimen


CAUTION!

Please set the microscope on stable table, desk, or other flat surface. Away from edge.

CARE & MAINTENANCE

1. Never leave the SV605 microscope close to fire, or hot places
2. When you are done using your microscope, remove any specimens left on the stage
3. Turn off the top and bottom illuminators when you are done using the microscope.
4. If you will not use the microscope for a long time (30 days or more), remove the batteries.
5. Always place the dust cap over the eyepiece when not in use or when being stored.
6. Store the microscope in a clean, dry place.
7. Clean the outside surfaces with a damp cloth.
8. Never clean optical surfaces with cloth or paper towels as they can scratch optical surfaces easily.
9. To remove fingerprints from optical surfaces, use a lens cleaning agent and lens tissue available at most photo outlets.
10. Never disassemble or clean internal optical surfaces. This should only be done by qualified technicians at the factory or other authorized repair facilities.

OVERVIEW

List and View of Accessory

SPECIFICATIONS

Model	SV605
Optical System	Finite Biological Microscope
Observation head	binocular head, 45 degree tilt, 360 degree rotation
Magnification	40X-1600X
Eyepiece	WF10X WF20X
Objective	4X 10X 40XS
Working Stage	90mmX90mm Square , Moving Scale , Disc Filters
Adjustment	Coarse Focuser & Fine Focuser (dual knobs) Focus range 15mm
Illuminator — Top	LED Lights
Illuminator — Bottom	LED Lights , Uses 2AA batteries (user supplied)
Accessories	Full Set (Accessory Chart)
Power supply	support 4*AA batteries and DC5V-1000MA input
Packing	Box
Dimension	18.5x14x29cm
Weight	1.75kg

SETTING UP & OPERATION

1. Take out the microscope and put the WF10X eyepiece into the eyepiece

2. Insert the power cord at the back or put the battery at the bottom to turn on the light source. Observe with eyepieces so that the field of view is fully illuminated. It can be adjusted by the light source adjustment knob.

3. Place the specimen slice on the working platform, clamp and adjust the position with a movable ruler, so that the observation object is in the center of the light-passing hole.

4. Lower the objective lens to the lowest position, and then slowly move it upward with coarse adjustment to obtain a relatively clear image, and then use fine adjustment to obtain a clearer image quality.

5. Turn the animal lens converter to adjust the different magnification objectives to complete the observation from low magnification to high magnification in sequence

6. The clarity of the observation object has a great relationship with the field of view illumination, and the upper and lower light sources must be adjusted. Rotate and adjust the color filter under the stage to change the field of view cable and color

Hand sliced production

Preparation of slices of animal and plant epidermis specimens

- For animals, take insect wings as an example: Use tweezers to tear a small piece of insect wings. Place it in the middle of the glass slide, and you can see the neatly arranged scales under the microscope. Wipe the scales directly on the glass slide, and you can observe the shape of a single scale. (The transparent wings can be observed directly) Apply concentrated salt water to the middle of the glass slide, and after drying, regular crystals can be observed under the microscope.
- For plants, take tomato pulp as an example: smear the pulp on the middle of the glass slide and drip with water. Then cover the cover glass, absorb the excess water with absorbent paper, and observe directly under the microscope.
- Some plant specimens, such as onion epidermis, root hairs, wheat awns, etc., can also be made into temporary specimens for observation by the above methods. However, some specimens, such as horizontal/vertical slices of plant stems, need to be peeled off by soaking, and then the above methods are used. Make it close to the slide to observe

NOTE! When cutting plant specimens, you need a slicer with a sharp blade inside. Be careful! When children make slices, they should be done under the guidance of adults!

- If glue is used instead of clean water, it is made by the above-mentioned method, and the cover fragment is covered, and the label is attached after it is dried. Then, a specimen section that is easy to observe and easy to preserve is made.
- The slide and cover fragment can be reused. The observation is complete. Wash the slides clean, dry them and store them.

Before using this device, read this guide which contains important operating instructions for safe usage and control for compliance with applicable standards and regulations.

FCC Requirements:

- Products authorized under Part 15 using SDoC or Certification require a label containing one of the following compliance statements
- (1)Receivers associated with licensed device service operations:
This device complies with part 15 of the FCC Rules. Operation is subject to the condition that this device does not cause harmful interference.
- (2)Stand-alone cable input selector switch:
This device complies with part 15 of the FCC Rules for use with cable television service.
- (3)All other devices:
•This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions:
(1) This device may not cause harmful interference, and
(2) this device must accept any interference received, including interference that may cause undesired operation.

CE Requirements:

•(Simple EU declaration of conformity) Hong Kong Svbonny

Technology Co.,Ltd declares that the equipment type is in compliance with the essential requirements and other relevant provisions of RED Directive 2014/30/EU and the ROHS Directive 2011/65/EU and the WEEE Directive 2012/19/EU; the full text of the EU declaration of conformity is available at the following internet address: www.svbony.com.

•Disposal

The crossed-out wheeled-bin symbol on your product, literature, or packaging reminds you that in the European Union, all electrical and electronic products, batteries, and accumulators (rechargeable batteries) must be taken to designated collection locations at the end of their working life. Do not dispose of these products as unsorted municipal waste. Dispose of them according to the laws in your area.

IC Requirements:

CAN ICES-3(B)/NMB-3(B)

Avoid Choking Hazard

Small Parts. Not for children under 3 years.

Approved Accessories

- This device meets the regulatory standards when used with the Svbonny accessories supplied or designated for the product.
- For a list of Svbonny-approved accessories for your item, visit the following website: <http://www.svbony.com>

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Warranty Card


Product Model:	
Purchasing Date:	
Defect Reason:	
Dealer Name:	
Telephone:	
User' s Name:	
User' s Address:	
User' s Email:	

Remarks:

1. This guarantee card should be kept by the user, no replacement if lost.
2. Most new products carry a one-year manufacturer' s warranty from the date of purchase.
3. The user can get warranty and after-sales service as below:
 - Contact the seller where you buy.
4. For warranty service, you will need to provide a receipt proof of purchase from the actual seller for verification

Exclusions from Warranty Coverage:

1. To any product damaged by accident.
2. In the event of misuse or abuse of the product or as a result of unauthorized alterations or repairs.
3. If the serial number has been altered, defaced, or removed.

CE FC RoHS

Hong Kong Svbonny Technology Co.,Ltd

Add: Unit B, 5th Floor, Gallo Commercial Building, 114-118
Lockhart Road, Wanchai, Hong Kong

Facebook: facebook.com/svbony

E-mail: info@svbony.com

Web: www.svbony.com

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